# Ultra compact pulsed Nd:YAG laser





#### **MAIN FEATURES**

- 1062, 532, 355, 266, 213 nm and 1.57 μm available
- Compact and portable, with quick umbilical disconnects
- Choice of resonators available to meet the need of demanding applications
- Integrated motorized variable attenuator and harmonics options
- Fiber coupling available at 1064 nm and 532 nm
- MIL-STD-810 tested to withstand harsh environments

## MAIN APPLICATIONS

- LIBS
- PUMPING
- FPD REPAIR
- LIDAR
- ABLATION
- PULSED LASER DEPOSITION
- PHOTOACOUSTIC IMAGING

#### Typical beam profiles



Ultra 50, Stable Near field @ 1064 nm



Ultra 50, GRM Far field @ 532 nm



Ultra 20, TEMOO Far field @ 355 nm



Ultra 50, Stable Near field @ 266 nm



Ultra 100, GRM Far field @ 213 nm

# www.quantel-laser.com

Many options and configurations are available. Please contact Lumibird to find the best match for your needs and compatibility between options.



Lumibird has locations across the globe that are available to provide support for any product, service or inquiry. Visit www.lumibird.com to connect with any of our global sites.



# **ULTRA**

# Ultra compact pulsed Nd:YAG laser

SPECIFICATI	ONS	ULTRA 20		ULTRA 50		ULTRA 100	
Resonator type		TEM00	Stable	Stable	GRM	Stable	GRM
Repetition rate (Hz)		1 to 20	1 to 50	1 to 20	20	1 to 20	20
Energy per pulse (mJ)	1064 nm	10	20	50		100	
	532 nm	6	12	30		55	
	355 nm	2	4	12		30	
	266 nm	1	4	10		25	15
	213 nm	-	-	-	-	4	On request
	1.57 µm	-	-	8	-	25	-
	1064 nm	<11	< 13	< 10	< 9	< 10.5	< 9
	532 nm	< 10	< 12	< 9	< 9	< 8	< 8.5
Pulse	355 nm	< 9	< 11	< 8	< 8	< 8	< 7.5
duration (ns) <sup>(1)</sup>	266 nm	< 9	< 11	< 8	< 8	< 8	< 8
	213 nm	-	-	-	-	< 7	-
	1.57 µm	-	-	< 10	-	< 9	-
Beam diameter (mm) (2)	1064 nm	1.3	2.5	3		4	
	1064 nm	< 2.5	< 6	< 7	< 1.5	< 8	< 1.5
	532 nm	< 1.5	< 5	< 6	< 1.5	< 7	< 1
Beam divergence (mrad) <sup>(3)</sup>	355 nm	< 1	< 4	< 5	< 1.2	< 5	< 1.5
	266 nm	< 1	< 4	< 7	< 1.5	< 4	< 1.5
	213 nm	-	-	-	-	< 3	-
	1.57 µm	-	-	< 12	-	< 12	-
Polarization ratio (%) <sup>(4)</sup>	1064 nm	> 99	> 99	> 98	> 98	> 98	> 98
Pulse to pulse energy stability (%) <sup>(5)</sup>	1064 nm	< 2	< 2.5	< 2	< 4	< 2	< 2
	532 nm	< 3.5	< 4	< 2.5	< 5	< 2.5	< 2.5
	355 nm	< 4	< 3	< 3	< 6	< 2	< 3
	266 nm	< 5	< 3	< 3	< 6	< 2	< 3
	213 nm	-	-	-	-	< 2	-
	1.57 µm	-	-	< 2	-	< 2	-



Power drift (%) (1)	1064 nm	± 5			
Pointing stability (µrad) <sup>(2</sup> )	All wavelengths	< 50			
Jitter (ns) (3)	All wavelengths	< ± 2			
Linewidth (cm-1) (4)	1064 nm	1			
(1) Over 8 hours, 18 °C < T < 28 °C, for $\Delta$ T < ± 5 °C					

(2) Full angle, 99 % of shots

(3) With respect to Q-Switch trigger

(4) Measured at FWHM

OTHER INF	ORMATION				
Power requirements		100-240 VAC, 50/60 Hz, single phase, 850 VA			
Cooling		Water to air			
Temperatures		Standard	Ethylene glycol option (EGW)		
	Operating	+ 10 °C to + 40 °C	+ 10 °C to + 40 °C		
	Storage	+ 5 °C to + 70 °C	- 30 °C to + 70 °C		
Laser head sealing		IP 66			
Vibration and shocks		Complies with MIL-STD-810			
Max. altitude (m)		3000			
Cable length (m)		3 (1)			
Flashlamps warranty		50 million shots <sup>(2)</sup>			
Weight (kg)	Laser head	0.9			
	Harmonic modules	0.5			
	Integrated cooling & electronics	14 (standard)	14. 5 (optional rack)		
1) Other lengths up to 15 m on request					

ther lengths up to 15 m on rec (2) 80 % of energy, or 1 year, whichever comes first

### Laser head 52 mm ULTRA 32 mm . 177 mm 76 mm **Remote Box**

04/21 REV B - Lumibird reserves the right to modify the specifications without prior notice.

 
 Image: 1.57 µm
 < 2</th>
 < 2</th>

 (1) Measured at FWHM with fast photodiode and 1 GHz oscilloscope
 (2) At the output of the laser
(3) Full angle, at 1/e<sup>2</sup> of the peak (4) Polarization is horizontal @ 1064 nm & 1.57 µm and vertical @ 532, 355, 266 & 213 nm (5) RMS, 99 % of shots **Integrated cooling OPTIONS** & electronics





Lumibird has locations across the globe that are available to provide support for any product, service or inquiry. Visit www.lumibird.com to connect with any of our global sites.

